

## 5 Function Electrical bed with battery back up



Low platform height 250 mm (10") from the floor  
• Safe working load 180 kg (396 lbs)  
• Electric backrest, knee break and height adjustment  
Battery back-up.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description -Low platform height 250 mm (10") from the floor  
• Safe working load 180 kg (396 lbs)  
• Electric backrest, knee break and height adjustment  
• Pro-Contour® Advance profiling feature has demonstrated low interface pressures which may help in the prevention of pressure injuries  
• Battery back-up.  
• Integrated mattress platform extension of 130 mm (5") with mattress retainer  
• Designed for use in private, residential and convalescent homes  
• Easy to disassemble into four parts for improved installation access  
Overall length 225 cm  
Overall width 105 cm

## BEDS: 5 Function Electrical bed with battery back up

---

Height adjustment 25 to 66 cm  
Mattress Platform length 202 cm  
Mattress Platform width 86 cm  
Backrest Tilt 69 degrees  
Footrest Tilt 32 degees  
Safe working load 180kg  
Maximum user weight 145kg