

## Advanced Chair Stair Evacuation Stretcher



Chair Stair stretcher is an emergency device utilised for evacuating injured or wheelchair bound people downstairs in emergency situations.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description - Manufactured from high strength aluminium alloy

- Utilised for evacuating patients downstairs in tall buildings
- Single or dual person operation
- Special wheel design provides excellent tracking and stability in confined spaces and in narrow stairwells
- Folds for easy storage
- Adjustable backrest ensures support for a range of patient heights
- Front and "free-loading" back handles allow for two-person operation
- Load bearing: 160kg