

## Trauma Stretcher



Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Description - Back rest:  $0-85^{\circ}\pm 5^{\circ}$ , controlled by gas spring.

- Knee break:  $0-30^{\circ}\pm 5^{\circ}$ , controlled by manual crank
- Lowest height:  $600\text{mm}\pm 10\text{mm}$
- Highest height:  $920\text{ mm}\pm 10\text{mm}$
- Tilt to back & forth which is controlled by the foot pedal on the two sides of the stretcher.
- Safety working loading: 2500N
- User environment : temperature:  $5^{\circ}\text{?}- 40^{\circ}\text{?}$
- The main frame of stretcher is made of galvanized steel tube coated with electrostatic powder spray.
- The top board is made of X-ray translucent material.